

Evaluation of LCMSMS Deuterium Scrambling in Clinically Significant Small Molecules

Authors: Joshua Cooper, Isil Dilek, and Uma Sreenivasan

Cerilliant Corporation, 811 Paloma Drive, Suite A, Round Rock, TX 78665

Introduction

Introduction and Objective:

LC-MS/MS is a powerful tool that brings numerous benefits to the clinical sample analysis arena. However, due to the complexity of the instrumentation there are some unique challenges that also accompany these benefits. Even following sample extraction and cleanup, matrix effects from the samples can cause interferences or impact ionization efficiency. Deuterium-labeled internal standards are the most common and prevalent labeled internal standards used to compensate for matrix effects. Some deuterium labeled compounds may exhibit hydrogen-deuterium scrambling/exchange in the collision cell which can impact MS/MS transition selection.

In this study we investigated numerous variables that potentially contribute to scrambling in order to ascertain reproducibility and impact on scrambling ratios: influences of different LC-MS systems (tandem quadrupole vs. quadrupole time-of-flight), matrix selection, concentration, with and without HPLC, collision energies, and deuterium placement in the internal standard. Numerous small molecules of clinical importance were investigated including: hydroxyvitamin D, testosterone, immunosuppressants, bath salts, and spice cannabinoids.

Materials and Procedures

LCMS Systems:



Waters Alliance UPLC-Xevo G2 Q-ToF
Agilent 1100 HPLC-6410 triple quad

Cerilliant Solution Standards Used:

25-Hydroxyvitamin D3, Cat# H-083
25-Hydroxyvitamin D3-d₆, Cat# H-074
25-Hydroxyvitamin D2, Cat# H-073
Testosterone, Cat# T-037
Testosterone-d₃, Cat# T-046
Testosterone-¹³C₃, Cat# T-037
(-)-Δ⁹-THC, Cat# T-005
(+)-Δ⁹-THC-d₃, Cat# T-003
(±)-11-Hydroxy-Δ⁹-THC, Cat# H-026
(±)-11-Hydroxy-Δ⁹-THC, Cat# H-041
(±)-11-nor-9-Carboxy-Δ⁹-THC, Cat# T-006
(±)-11-nor-9-Carboxy-Δ⁹-THC-d₃, Cat# T-004
(±)-11-nor-9-Carboxy-Δ⁹-THC-d₆, Cat# T-007
Cannabinol, Cat# C-045
Cannabinol-d₃, Cat# C-084
JWH-018 4-Hydroxypentyl metabolite, Cat# S-035
JWH-018 4-Hydroxypentyl metabolite-d₅, Cat# S-039
JWH-073 3-Hydroxybutyl metabolite, Cat# S-037
JWH-073 3-Hydroxybutyl metabolite-d₅, Cat# S-040
3,4-MDPV HCl, Cat# M-146
3,4-MDPV-d₃ HCl, Cat# M-146
Ethylone HCl, Cat# E-071
Ethylone-d₅ HCl, Cat# E-072
Butylone HCl, Cat# B-045
Butylone-d₃ HCl, Cat# B-046
Mephedrone HCl, Cat# M-138
Mephedrone-d₃ HCl, Cat# M-139
Methylone HCl, Cat# M-140
Methylone-d₃ HCl, Cat# M-140
Everolimus-d₄, Cat# E-070
Mycophenolic acid, Cat# M-106
Mycophenolic acid-d₃, Cat# M-137

Serum Extraction:
200µl of sample in serum + 200µl of methanol, vortexed to mix.
Added 1 ml of heptane, vortexed for 30sec,
Centrifuged for 4min at 3000rpm
900µl of top layer dried under nitrogen
Reconstituted in 100µl of ethanol

Comparisons of 25-Hydroxyvitamin D Deuterium Scrambling



Labeled 25-Hydroxyvitamin D2 and D3 Scrambling in Serum

| Compound | Label | System | Concentration µg/mL | Transition d _{n-1} | Transition d _n | Scrambling % d _{n-1} /d _n |
|----------------------|----------------|---------|---------------------|-----------------------------|---------------------------|---|
| 25-Hydroxyvitamin D2 | d ₃ | Xevo G2 | 2 | 398→379 | 398→380 | 28.6 |
| | | | 0.2 | 398→379 | 398→380 | 35.4 |
| | | 6410 | 5 | 416→397 | 416→398 | 2.8 |
| | | | | 416→379 | 416→380 | 19.7 |
| | | | 50 | 398→379 | 398→380 | 30.4 |
| | | | | 416→397 | 416→398 | 2.8 |
| | d ₆ | 6410 | 5 | 419→400 | 419→401 | 2 |
| | | | | 419→382 | 419→383 | 8.8 |
| | | 50 | 401→382 | 401→383 | 5.9 | |
| | | | 419→400 | 419→401 | 2 | |
| | | | 419→382 | 419→383 | 9 | |
| | | | 401→382 | 401→383 | 5.4 | |
| 25-Hydroxyvitamin D3 | d ₆ | 6410 | 2.5 | 407→388 | 407→389 | 4 |
| | | | 407→370 | 407→371 | 18.8 | |
| | | | | 389→370 | 389→371 | 9.2 |

Transitions Comparisons for Native and Labeled 25-Hydroxyvitamin D2 and D3 in EtOH on 6410

| Parent → Water loss | | | | | |
|----------------------|----------------|---------------------|-----------------------------|---------------------------|---|
| Compound | Label | Concentration µg/mL | Transition d _{n-1} | Transition d _n | Scrambling % d _{n-1} /d _n |
| 25-Hydroxyvitamin D2 | d ₃ | 100 | 416→397 | 416→398 | 2.9 |
| | d ₆ | 100 | 419→400 | 419→401 | 2 |
| | native | 50 | 413→394 | 413→395 | 0.5 |
| 25-Hydroxyvitamin D3 | d ₆ | 50 | 407→388 | 407→389 | 4 |
| | native | 100 | 401→382 | 401→383 | 0.5 |


| Parent → 2 Water losses | | | | | |
|-------------------------|----------------|---------------------|-----------------------------|---------------------------|---|
| Compound | Label | Concentration µg/mL | Transition d _{n-1} | Transition d _n | Scrambling % d _{n-1} /d _n |
| 25-Hydroxyvitamin D2 | d ₃ | 100 | 416→379 | 416→380 | 19.5 |
| | d ₆ | 100 | 419→382 | 419→383 | 8.9 |
| | native | 50 | 413→376 | 413→377 | 0.5 |
| 25-Hydroxyvitamin D3 | d ₆ | 50 | 407→370 | 407→371 | 18.9 |
| | native | 100 | 401→364 | 401→365 | 0.3 |

| Water Loss → 2 Water losses | | | | | |
|-----------------------------|----------------|---------------------|-----------------------------|---------------------------|---|
| Compound | Label | Concentration µg/mL | Transition d _{n-1} | Transition d _n | Scrambling % d _{n-1} /d _n |
| 25-Hydroxyvitamin D2 | d ₃ | 100 | 398→379 | 398→380 | 30.4 |
| | d ₆ | 100 | 401→382 | 401→383 | 5.4 |
| | native | 50 | 398→376 | 398→377 | 0.4 |
| 25-Hydroxyvitamin D3 | d ₆ | 50 | 389→370 | 389→371 | 11.2 |
| | native | 100 | 383→364 | 383→365 | 0.3 |

Notes: 25-Hydroxy D2-d₆ water loss→2 water loss has same transition as 25-Hydroxyvitamin D3 parent→water loss. Can be problem if compounds are not well resolved chromatographically.

Selection of Transitions Greatly Impacts Observed Scrambling

Investigation of Testosterone Scrambling




Major transitions are:
Native: 289→97 & 289→109
Testosterone-d₃: 292→97 & 292→109
Testosterone-¹³C₃: 292→100 & 292→112
No scrambling at major transitions

Testosterone Scrambling at Minor Transitions

| Label | Method | Instrument | Concentration µg/ml | Transitions D _{n-1} or ¹³ C _{n-1} | Transitions D _n or ¹³ C _n | *Scrambling % D _{n-1} / D _n |
|----------------|----------|------------|---------------------|--|--|---|
| d ₃ | Infusion | Q-ToF | 10 | 292→255 | 292→256 | 31.9 |
| | | | 100 | | | 36.5 |
| | | | 10 | | | 35.7 |
| | LC | 6410 | 10 | 292→255 | 292→256 | 37.7 |
| | | | 100 | | | 36.3 |
| | | | 10 | | | 0.1 |
| native | | | 100 | 289→252 | 289→253 | 0.0 |


* or Scrambling % ¹³C_{n-1} / ¹³C_n

Testosterone Chromatograms on Xevo G2




Investigation of Scrambling in Spice Cannabinoids, Bath Salts, and Immunosuppressants


(±)-11-nor-9-Carboxy-Δ⁹-THC Scrambling at m/z 327



JWH-073 3-Hydroxybutyl Metabolite Scrambling at m/z 216



Everolimus Scrambling at m/z 389



Spice Cannabinoids, Bath Salts, and Immunosuppressants Scrambling Comparison using Xevo G2

| Compound | Label | Polarity | Collision Energy | Transition(s) d _n | Scrambling % d _{n-1} /d _n |
|--|----------------|----------|------------------|------------------------------|---|
| (-)-Δ ⁹ -THC | d ₃ | pos | 25 | 318→262, 196 | 0 |
| | native | pos | 25 | 315→259, 193 | 0 |
| (+)-Δ ⁹ -THC | d ₃ | neg | 30 | 318→262, 196 | can't determine |
| | native | neg | 30 | 315→259, 193 | can't determine |
| (±)-11-Hydroxy-Δ ⁹ -THC | d ₃ | pos | 15 | 334→316 | 3.12 |
| | native | pos | 15 | 331→313 | 0 |
| (±)-11-Hydroxy-Δ ⁹ -THC | d ₃ | pos | 25 | 334→any | can't determine |
| | native | pos | 25 | 331→any | can't determine |
| (±)-11-nor-9-Carboxy-Δ ⁹ -THC | d ₆ | neg | 30 | 352→254 | 0 |
| | d ₃ | neg | 30 | 346→248 | 0 |
| | native | neg | 30 | 343→245 | 0 |
| (±)-11-nor-9-Carboxy-Δ ⁹ -THC | d ₆ | neg | 20 | 352→334 | 0 |
| | d ₃ | neg | 20 | 346→328 | 0 |
| | native | neg | 20 | 343→325 | 0 |
| (±)-11-nor-9-Carboxy-Δ ⁹ -THC | d ₆ | pos | 20 | 354→336 | 51.91 |
| | d ₃ | pos | 20 | 348→330 | 0 |
| | native | pos | 20 | 345→327 | 0 |
| (±)-11-nor-9-Carboxy-Δ ⁹ -THC | d ₆ | pos | 20 | 354→308 | 48.88 |
| | d ₃ | pos | 20 | 348→302 | 0 |
| | native | pos | 20 | 345→299 | 0 |
| Cannabinol | d ₃ | neg | 30 | 316→248 | 0 |
| | native | neg | 30 | 313→245 | 0 |
| Cannabinol | d ₃ | pos | 20 | 318→262, 196 | 0 |
| | native | pos | 20 | 315→259, 193 | 0 |
| JWH-018 4-Hydroxypentyl metabolite | d ₅ | pos | 20 | 363→345 | 22.83 |
| | native | pos | 20 | 358→340 | 0 |
| JWH-018 4-Hydroxypentyl metabolite | d ₅ | pos | 20 | 363→155 | 0 |
| | native | pos | 20 | 358→155 | 0 |
| JWH-073 3-Hydroxybutyl metabolite | d ₅ | pos | 20 | 349→221 | 25.44 |
| | native | pos | 20 | 344→216 | 0 |
| JWH-073 3-Hydroxybutyl metabolite | d ₅ | pos | 20 | 363→155 | 0 |
| | native | pos | 20 | 358→155 | 0 |
| 3,4-MDPV HCl | d ₈ | pos | 15 | 284→134 | 0 |
| | native | pos | 15 | 284→126 | 0 |
| Ethylone HCl | d ₅ | pos | 15 | 227→209 | 0 |
| | native | pos | 15 | 222→204 | 0 |
| Butylone HCl | d ₃ | pos | 15 | 225→209, etc | 0 |
| | native | pos | 15 | 222→204, etc | 0 |
| Mephedrone HCl | d ₃ | pos | 10 | 181→163 | 0 |
| | native | pos | 10 | 178→160 | 0 |
| Methylone HCl | d ₃ | pos | 10 | 211→163 | 0 |
| | native | pos | 10 | 208→160 | 0 |
| Methylone HCl | d ₃ | pos | 10 | 211→135 | 0 |
| | native | pos | 10 | 208→132 | 0 |
| Everolimus | d ₄ | pos | 60 | 984→393 | 1.80 |
| | native | pos | 60 | 980→389 | 0 |
| Mycophenolic acid | d ₃ | neg | 15 | 322→278 | 0 |
| | native | neg | 15 | 319→275 | 0 |

CONCLUSIONS

- Scrambling was observed for several of the analytes at select transitions. In all cases, scrambling was mitigated or eliminated by optimizing instrument conditions and transition selection.
- Awareness of potential scrambling is important for proper internal standard selection.
- Scrambling was observed on both the Agilent 6410 triple quadrupole and the Waters Xevo G2 Q-ToF to approximately the same degree. For a specific transition, scrambling ratios were consistent between solvent and serum. No matrix effects on scrambling.
- Direct infusion can provide rapid and accurate determination of scrambling ratios. Infusion and chromatographic injection results were consistent.
- Scrambling may be mitigated or eliminated by altering instrument conditions and transition selection. However, there is a need to consider potential impact of scrambling on transitions chosen for optimal sensitivity.
- Deuterium-labeled internal standards are a viable option for LC-MS/MS analysis with selection of the appropriate transition. Deuterated standards can be more cost effective than ¹³C labeled internal standards, more widely available and with lower cost per test. ¹³C labeled internal standards are most effective when deuterium scrambling issues can not be resolved.