

Measurement of Thyroxine (T4), Triiodothyronine (T3) and Reverse Triiodothyronine (rT3) by Liquid Chromatography with Online Sample Cleanup-Tandem Mass Spectrometry in Negative Mode

Dorothy Yang, Anabel Fandino, Agilent Technologies, Inc. Santa Clara, CA

MSACL 2013
Poster # 4



Introduction

Liquid chromatography and triple quadrupole mass spectrometry (LC-MS/MS) has become an essential tool for small molecules quantitation due to its high sensitivity and specificity, excellent reproducibility and the ability to perform simultaneous analysis of multiple analytes. The accurate and precise measurement of thyroid hormones in blood is imperative in thyroid function monitoring. Quantitation of thyroid hormones in serum or plasma is challenging due to their low levels under normal physiological conditions. We have developed a sensitive and reliable LC-MS/MS method on Agilent 6490 QQQ that utilizes JetStream ionization source and dual ion funnel technologies. These innovative technologies enable more efficient ion generation and ion sampling that result in the sensitive detection of Thyroxine (T4), Triiodothyronine (T3) and Reverse Triiodothyronine (rT3) in serum. We also explore the online sample cleanup to minimize the workload of sample preparation

