



**Cerilliant**<sup>®</sup>

Analytical Reference Standards

a SIGMA-ALDRICH<sup>®</sup> company

## Recovery of Codeine-6-Glucuronide with Beta-Glucuronidase

Several customers have reported lower than expected recovery of Codeine-6-glucuronide with beta-glucuronidase.

Some of our customers have resolved the issue and provided suggestions. The problem appears related to substrate specific enzyme reaction rate and/or capacity, with better reaction rate and conversion for morphine-3-glucuronide than opiate-6-glucuronides. They have observed better results using a 4 fold excess of enzyme relative to glucuronide over their normal procedure. This can be achieved either by reducing the concentration/volume of the glucuronide or increasing the concentration/volume of the enzyme. For LC/MS/MS applications lowering the concentration of glucuronide, buffer and IS to 1/4th was recommended due to lower sensitivity of the LC/MS/MS technique relative to GC/MS. Stereochemical position of the glucuronide or activity or type of glucuronidase could be factors that affect reaction rate and recovery.

811 PALOMA DRIVE, SUITE A, ROUND ROCK, TEXAS 78665

PHONE 512/238-9974 | 800/848-7837 | FAX 512/238-9129 | 800/654-1458 | [www.cerilliant.com](http://www.cerilliant.com)